InnoVet Select HF

HIGH FREQUENCY RADIOGRAPHIC SYSTEM



"THE MOST TRUSTED NAME IN VETERINARY RADIOGRAPHY"

Designed by veterinarians for veterinarians, the InnoVet Select HF has "All the little things in mind" and delivers uncompromising functionality.

INNOVET MEANS INTEGRITY

InnoVet revolutionized veterinary radiography in 1991 with the first x-ray system designed by veterinarians. Since then, InnoVet has manufactured more than 6,000 veterinary x-ray systems, making it the #1 veterinary x-ray manufacturer in the

United States.

Add to that, the industry's best warranty and it's no wonder InnoVet is the most trusted name in veterinary radiography.



LOCAL SUPPORT

Capital equipment is a long-term investment that deserves personal attention from a service professional who understands your practice.

Regional InnoVet dealers are factory-trained to install and service your x-ray systems with immediate on site and remote support. Our trusted support network means image quality will never slow down your practice.

INDUSTRY BEST WARRANTY

Our 7-year limited parts warranty secures your investment today and guarantees your investment value tomorrow.

DIGITAL READY - RIGHT FOR TODAY OR TOMORROW

Our systems are compatible with most digital imaging systems.

2901 W. Lawrence Ave. • Chicago, IL 60625 800-972-9776 • 773-588-2444 Fax 773-588-3424

www.innovet4vets.com



System Specifications

HIGH FREQUENCY GENERATOR

- 20 KHz High Frequency technology
- 30 kW rating
- Two-point technique selection (kVp, mAs)
- kVp range: 40 125, in 1 kV increments
- mA stations: 50s, 100s, 100L, 300L
- Digital kVp, mA, and mAs readout
- Tilt-out control panel
- Built-in tube protection and self diagnostics
- Adjustable volume, audible exposure signal
- Automatic generator "Shut Off" circuit
- Tilt-out control panel
- UL classified

TABLE / TUBESTAND

- Enclosed, integrated radiographic table
- Stationary table top: 53 in. x 24 in.
- Welded base construction
- · Transformer mounting shelf
- · Grid cabinet with Simplex cassette tray
- Grid cabinet / tubestand interlock
- Grid 8:1 ratio, 103 LPI
- Animal restraining device
- Angulation dial and deluxe operator handgrips
- Tubestand handles

X-RAY TUBF

• 1.0-2.0 mm F.S., 140,000 H.U. with 90° horns

• 12 ft. high voltage cables with federal terminations

COLLIMATOR

- Manual collimator with light field and swivel mount
- Field light "on" time selectable 30, 60, 90 or 120 seconds
- Radiation leakage well below Federal standards for human medical systems

System Options

- Anatomical programmed technique selection
- 500 mA generator upgrade
- Two-way float top table, 58 in. x 32 in.
- Four-way float top table, 67 in. x 32 in.
- Table top fluid trough
- "Hands-free" collimator light system