




ZEN-PX4 X-RAY



The ZEN-PX4 is the newest addition to our dental X-ray offering. Handheld X-rays are becoming a popular choice for veterinarians looking to expand dentistry in their practice. The ZEN-PX4 features a small and lightweight design making it easy to use and store when not in use. Building off the success of the PX2 handheld X-ray, the PX4 is equipped with an upgraded battery rated for up to 400 exposures on one charge.

Handheld X-rays are a great choice for practices both new to and familiar with dental imaging. Complete mobility allows these units to be effortlessly deployed to any treatment table in your veterinary practice as dentistry continues to grow. The next generation of handheld X-rays begins with the ZEN-PX4.

Specifications:

Weight: 3.3lbs

Battery: 22.2V Li-Polymer Battery
(400 Exposures @ 0.5 Seconds)

Dimensions: 6"x7"x10"